

# Steps Involved In Cleaning A NEBULIZER



Place the mask or mouthpiece, top piece, and medicine cup on the dishwasher's top shelf or hand-wash them in warm soapy water.



Let the units air-dry in a cool and dry place.

Once a week, deep-cleanse these segments by soaking them in a white vinegar-water solution for 30 minutes or as recommended in the user's manual.



Never wash the tubing and the compressor. Clean them with a disinfectant wipe.

Replace the compressor's air filter every six months or as the manufacturer recommends.



Clean the nebulizer after every use.

#### Reference

How To Use A Nebulizer; American Lung Association